

Kollidon® SR

PRD-No.: 30071321

RegXcellence®

Product Specification

Workflow No.: DAWF-2024-0299

Revision 8

Effective from: 01 Jul 2024

Status: final

Test parameter	Requirements	Test method
Appearance	white to slightly yellowish, free flowing powder	visual
Identification		
Identification (IR)	Must comply	PM/00511
Identification (sulfate)	Must comply	Ph.Eur. ("sodium lauril sulfate", ID test C)
Purity tests		
Content of Povidone (calculated: Nitrogen / 0.126)	18.0 to 21.0 g/100g (%)	PM/02047
Content of Polyvinylacetate (calculated: Saponification value x 0.1534)	76.0 to 82.0 g/100g (%)	PM00529QC
Content of sodium lauril sulfate	≤ 0.7 g/100g (%)	PM02385QC
Content of Silica (SiO ₂) (calculated: Si x 2.142)	0.4 to 0.7 g/100g (%)	PM/02049
Sulphated ash / Residue on ignition (weight sample: 10 - 20 g)	≤ 1.7 g/100g (%)	Ph.Eur., 2.4.14 (equivalent to USP-NF <281> and JP 2.44)
Residual solvents		
Acetic Acid, Acetate (class 3)	≤ 0.5 g/100g (%)	PM01926QC
pH value (suspension 10 % in water)	4.0 to 5.0	USP-NF, <791>
Vinyl acetate	≤ 100 mg/kg (ppm)	PM02232QC
Loss on drying	≤ 4.0 g/100g (%)	PM00993QC
Other tests		
Microbiological enumeration test		Ph.Eur. 2.6.12
Total aerobic microbial count(TAMC)	≤ 100 CFU/g	
Total combined yeasts and moulds count (TYMC)	≤ 10 CFU/g	
Test for specified microorganisms		Ph.Eur. 2.6.13
Escherichia coli	Absent/g	
Pseudomonas aeruginosa	Absent/g	
Staphylococcus aureus	Absent/g	
Bile-tolerant gram-negative bacteria	≤ 10 CFU/g	

For qualification, not for submission purposes. Submission-relevant specifications are compendia-specific, part of our regulatory documentation, and can be requested via [RegXcellence \(https://virtualpharmaassistants.basf.com\)](https://virtualpharmaassistants.basf.com)

BASF'S DISCLAIMERS AND WAIVERS ARE SET FORTH ON THE LAST PAGE.

If you have any further questions or need additional support, please contact your BASF sales representative.

Kollidon® SR

PRD-No.: 30071321

Product Specification

RegXcellence®

Workflow No.: DAWF-2024-0299

Revision 8

Test parameter

Requirements

Test method

The product meets the requirements of this specification.

This is a computer generated document. It is valid without signature.

Kollidon® SR

PRD-No.: 30071321

Product Specification

RegXcellence®

Workflow No.: DAWF-2024-0299

Revision 8

Except as expressly agreed otherwise in writing in a supply contract or other written agreement between you and BASF, this document and any information provided herein, does not constitute a legally binding obligation of BASF and has been prepared in good faith and is believed to be accurate as of the date of issuance. BASF rejects any obligation to, and will not, automatically update this document and any information provided herein, unless required by applicable law.

Unless agreed otherwise in writing, you are responsible for confirming that you have retrieved the most current version of this document from BASF. Any information provided herein can be changed at BASF's sole discretion anytime and neither this document nor the information provided herein may be relied upon to satisfy any obligations you may have to undertake your own inspections and evaluations.

To the fullest extent not prohibited by the applicable laws, BASF EXPRESSLY DISCLAIMS ALL REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES OF ANY KIND, WHETHER EXPRESS OR IMPLIED, WRITTEN OR ORAL, BY FACT OR LAW, INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTIES, REPRESENTATIONS OR CONDITIONS OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, SATISFACTORY QUALITY, NON-INFRINGEMENT, AND ANY REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES, ARISING FROM STATUTE, COURSE OF DEALING OR USAGE OF TRADE AND BASF HEREBY EXPRESSLY EXCLUDES AND DISCLAIMS ANY LIABILITY RESULTING FROM OR IN CONNECTION WITH THIS DOCUMENT OR ANY INFORMATION PROVIDED HEREIN, INCLUDING, WITHOUT LIMITATION, ANY LIABILITY FOR ANY DIRECT, CONSEQUENTIAL, SPECIAL, OR PUNITIVE DAMAGES RELATING TO OR ARISING THEREFROM, EXCEPT IN CASES OF

- (I) DEATH OR PERSONAL INJURY,
- (II) BASF'S OR ITS AGENTS AND ASSISTANTS WILLFUL MISCONDUCT, FRAUD OR FRAUDULENT MISREPRESENTATION OR
- (III) ANY MATTER IN RESPECT OF WHICH IT WOULD BE UNLAWFUL FOR BASF TO EXCLUDE OR RESTRICT LIABILITY UNDER THE APPLICABLE LAW.

If you have any further questions or need additional support, please contact your BASF sales representative

Attention! Please note that the certificates of analysis are also conveniently available on your BASF online portal.

PFIZER INC

2025-12-11

Certificate No 48

Page 1 of 2

445 EASTERN POINT RD, Bldg 98

EN - Nutrition & Health

06340-5157 GROTON

ENT/QA

QAL-Release@basf.com

USA

Inspection Certificate 3.1 according to EN 10204

Kollidon® SR

Material

52622254

Order

0164269538 000010

Delivery

0149550935 000010

Batch

12657788Q0

Batch/Qty

20.000 KG

Total

20.000 KG

20KG Plastic drums

Purchase order/Customer material

9501093114 / 1400

NO CODE NUMBER

Manufacturing date : 06/17/2025

Release date : 08/29/2025

Retest date : 06/16/2029

Characteristic	Unit	Value	Lower Limit	Upper Limit
Appearance		Conforms	white to slightly yellowish, free flowing powder	
Identification (IR)		Conforms	must conform	
Identification (sulfate)		Conforms	must conform	
Content of povidone	g/100 g	19.3	18.0	21.0
Content of polyvinylacetate	g/100 g	77.6	76.0	82.0
Content of sodium lauril sulfate	g/100 g	0.7		0.7
Content of silica (SiO ₂)	g/100 g	0.6	0.4	0.7
Sulphated ash / Residue on ignition	g/100 g	1.4		1.7
Residual solvent Acetic Acid, acetate	g/100 g	0.3		0.5
pH-value (suspension 10 % in water)		4.7	4.0	5.0
Vinyl acetate	mg/kg	< 30		100
Loss on drying	g/100 g	2.3		4.0
Total aerobic microbial count (TAMC)	CFU/g	< 10		100
Total combined yeasts/moulds count (TYMC)	CFU/g	< 10		10
Escherichia coli	CFU/g	0		0
Pseudomonas aeruginosa	CFU/g	0		0
Staphylococcus aureus	CFU/g	0		0
Bile-tolerant gram-negative bacteria	CFU/g	< 1		10

Manufacturer: BASF SE

Carl-Bosch-Str. 38

67056 Ludwigshafen

Germany

The above-mentioned data describes the relevant parameters of the product at the time of issuing this certificate of analysis (CoA). The data is controlled at regular intervals as part of our quality assurance program and is provided for reference only. The contractually agreed quality of the product at the time of passing of risk is exclusively determined by our product specification. Unless specifically agreed by the Parties in writing signed by authorized representatives, no suitability for a particular use shall be presumed or implied based on the CoA and the customer is solely responsible for determining and verifying the suitability of the use of the product for any purpose.

Attention! Please note that the certificates of analysis are also conveniently available on your BASF online portal.

PFIZER INC

2025-12-11

Certificate No 48

Page 2 of 2

445 EASTERN POINT RD, Bldg 98

EN - Nutrition & Health

06340-5157 GROTON

ENT/QA

QAL-Release@basf.com

USA

Inspection Certificate 3.1 according to EN 10204

Kollidon® SR

Material

52622254

Order

0164269538 000010

Delivery

0149550935 000010

Batch

12657788Q0

Batch/Qty

20.000 KG

Total

20.000 KG

20KG Plastic drums

Purchase order/Customer material

9501093114 / 1400

NO CODE NUMBER

Characteristic	Unit	Value	Lower Limit	Upper Limit
----------------	------	-------	-------------	-------------

Co-Manufacturer: Schwarzwaldmilch GmbH
Offenburg
Germany

BASF SE

Quality Assurance - Release

The above-mentioned data describes the relevant parameters of the product at the time of issuing this certificate of analysis (CoA). The data is controlled at regular intervals as part of our quality assurance program and is provided for reference only. The contractually agreed quality of the product at the time of passing of risk is exclusively determined by our product specification. Unless specifically agreed by the Parties in writing signed by authorized representatives, no suitability for a particular use shall be presumed or implied based on the CoA and the customer is solely responsible for determining and verifying the suitability of the use of the product for any purpose.

Kollidon® SR

PRD-No.: 30071321

Compliance with International Standards

RegXcellence®

Compliance with International Standards

- ☐ GMP for APIs (ICH Q7 Guideline)
- ☒ GMP for Excipients (Joint IPEC-PQG GMP Guide for Pharmaceutical Excipients/ EXCiPACT)
- ☐ HACCP
- ☐ ISO 9001
- ☒ co-manufacturer Quality status: GMP for Excipients
- ☐ Others:

Compendial compliance

- ☐ European Pharmacopoeia (Ph. Eur.)
- ☐ US Pharmacopoeia (USP-NF)
- ☐ Japanese Pharmacopoeia (JP)
- ☐ Japanese Pharmaceutical Excipients (JPE)
- ☐ Chinese Pharmacopoeia (ChP)
- ☐ Indian Pharmacopoeia (IP)
- ☐ Others:

For qualification, not for submission purposes – further compliance documentation can be requested via [RegXcellence® \(https://virtualpharmaassistants.basf.com\)](https://virtualpharmaassistants.basf.com)

Kollidon® SR

PRD-No.: 30071321

Compliance with International Standards

RegXcellence®

This document, or any information provided herein does not constitute a legally binding obligation of BASF and has been prepared in good faith and is believed to be accurate as of the date of issuance. Unless expressly agreed otherwise in writing in a supply contract or other written agreement between you and BASF:

(a) To the fullest extent not prohibited by the applicable laws, BASF EXPRESSLY DISCLAIMS ALL OTHER REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES OF ANY KIND, WHETHER EXPRESS OR IMPLIED, WRITTEN OR ORAL, BY FACT OR LAW, INCLUDING ANY IMPLIED WARRANTIES, REPRESENTATIONS OR CONDITIONS OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, SATISFACTORY QUALITY, NON-INFRINGEMENT, AND ANY REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES, ARISING FROM STATUTE, COURSE OF DEALING OR USAGE OF TRADE and BASF HEREBY EXPRESSLY EXCLUDES AND DISCLAIMS ANY LIABILITY RESULTING FROM OR IN CONNECTION WITH THIS DOCUMENT OR ANY INFORMATION PROVIDED HEREIN, including, without limitation, any liability for any direct, consequential, special, or punitive damages relating to or arising therefrom, except in cases of:

- (i) death or personal injury to the extent caused by BASF's sole negligence;
- (ii) BASF's willful misconduct, fraud or fraudulent misrepresentation; or
- (iii) any matter in respect of which it would be unlawful for BASF to exclude or restrict liability under the applicable laws;

(b) Any information provided herein can be changed at BASF's sole discretion anytime and neither this document nor the information provided herein may be relied upon to satisfy from any and all obligations you may have to undertake your own inspections and evaluations;

(c) BASF rejects any obligation to, and will not, automatically update this document and any information provided herein, unless required by applicable law; and

(d) The user is responsible for confirming that the user has retrieved the most current version of this document from BASF as appropriate.

Pharmacopoeia Compliance Statement

To whom it may concern

In line with the General Notices of Ph.Eur., JP, and USP, it is not required to test each batch for each monograph requirement of the applicable monograph, or to apply solely compendial methods for product release. The FDA Compliance Policy Guide* accepts such compendial standards.

Instead, the assurance of pharmacopoeia quality may be achieved by alternative control strategies. For example, process validation and appropriate in-process controls can justify the omission of tests listed in a monograph, including identity testing of any kind. Appropriate method validation allows for a reduced test program and compliance statements with different pharmacopoeias at the same time.

Therefore, BASF may apply reduced testing and validated alternative test methods for product release.

Upcoming changes to applicable, compendial monographs are monitored and will be implemented as they become official. Compliance with any cited, current monograph at the time of product release is stated on the Certificate of Analysis.

* U.S. Food and Drug Administration, CPG Sec. 420.400, "Performance of Tests for Compendial Requirements on Compendial Products" (Oct. 1, 1980).

This is a computer-generated document. It is valid without signature.

BASF'S DISCLAIMERS AND WAIVERS ARE SET FORTH ON THE LAST PAGE

If you have any further questions or need additional support, please contact your BASF sales representative.

Pharmacopoeia Compliance Statement

RegXcellence®

This document, or any information provided herein does not constitute a legally binding obligation of BASF and has been prepared in good faith and is believed to be accurate as of the date of issuance. Unless expressly agreed otherwise in writing in a supply contract or other written agreement between you and BASF:

- (a) To the fullest extent not prohibited by the applicable laws, BASF EXPRESSLY DISCLAIMS ALL OTHER REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES OF ANY KIND, WHETHER EXPRESS OR IMPLIED, WRITTEN OR ORAL, BY FACT OR LAW, INCLUDING ANY IMPLIED WARRANTIES, REPRESENTATIONS OR CONDITIONS OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, SATISFACTORY QUALITY, NON-INFRINGEMENT, AND ANY REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES, ARISING FROM STATUTE, COURSE OF DEALING OR USAGE OF TRADE and BASF HEREBY EXPRESSLY EXCLUDES AND DISCLAIMS ANY LIABILITY RESULTING FROM OR IN CONNECTION WITH THIS DOCUMENT OR ANY INFORMATION PROVIDED HEREIN, including, without limitation, any liability for any direct, consequential, special, or punitive damages relating to or arising therefrom, except in cases of
 - (i) death or personal injury to the extent caused by BASF's sole negligence,
 - (ii) BASF's willful misconduct, fraud or fraudulent misrepresentation or
 - (iii) any matter in respect of which it would be unlawful for BASF to exclude or restrict liability under the applicable laws;
- (b) Any information provided herein can be changed at BASF's sole discretion anytime and neither this document nor the information provided herein may be relied upon to satisfy from any and all obligations you may have to undertake your own inspections and evaluations;
- (c) BASF rejects any obligation to, and will not, automatically update this document and any information provided herein, unless required by applicable law; and
- (d) The user is responsible for confirming that the user has retrieved the most current version of this document from BASF as appropriate

Kollidon® SR
PRD-No.: 30071321

Statement: BSE-TSE (Bovine Spongiform Encephalopathy - Transmissible Spongiform Encephalopathy)

The above mentioned product is manufactured without any material of animal or human origin. The product meets the requirements of the „Note for Guidance on Minimising the Risk of Transmitting Animal Spongiform Encephalopathy Agents via Human and Veterinary Medicinal Products” (EMA/410/01 current version).

This document has been created digitally and is valid without a signature.

BASF'S DISCLAIMERS AND WAIVERS ARE SET FORTH ON THE LAST PAGE

If you have any further questions or need additional support, please contact your BASF sales representative.

Kollidon® SR

PRD-No.: 30071321

Statement: BSE-TSE (Bovine Spongiform Encephalopathy - Transmissible Spongiform Encephalopathy)

This document, or any information provided herein does not constitute a legally binding obligation of BASF and has been prepared in good faith and is believed to be accurate as of the date of issuance. Unless expressly agreed otherwise in writing in a supply contract or other written agreement between you and BASF:

(a) To the fullest extent not prohibited by the applicable laws, BASF EXPRESSLY DISCLAIMS ALL OTHER REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES OF ANY KIND, WHETHER EXPRESS OR IMPLIED, WRITTEN OR ORAL, BY FACT OR LAW, INCLUDING ANY IMPLIED WARRANTIES, REPRESENTATIONS OR CONDITIONS OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, SATISFACTORY QUALITY, NON-INFRINGEMENT, AND ANY REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES, ARISING FROM STATUTE, COURSE OF DEALING OR USAGE OF TRADE and BASF HEREBY EXPRESSLY EXCLUDES AND DISCLAIMS ANY LIABILITY RESULTING FROM OR IN CONNECTION WITH THIS DOCUMENT OR ANY INFORMATION PROVIDED HEREIN, including, without limitation, any liability for any direct, consequential, special, or punitive damages relating to or arising therefrom, except in cases of:

- (i) death or personal injury to the extent caused by BASF's sole negligence;
- (ii) BASF's willful misconduct, fraud or fraudulent misrepresentation; or
- (iii) any matter in respect of which it would be unlawful for BASF to exclude or restrict liability under the applicable laws;

(b) Any information provided herein can be changed at BASF's sole discretion anytime and neither this document nor the information provided herein may be relied upon to satisfy from any and all obligations you may have to undertake your own inspections and evaluations;

(c) BASF rejects any obligation to, and will not, automatically update this document and any information provided herein, unless required by applicable law; and

(d) The user is responsible for confirming that the user has retrieved the most current version of this document from BASF as appropriate

Kollidon® SR
PRD-No.: 30071321

Statement: Melamine free

According to our knowledge of the production process, raw materials, and equipment used, the above mentioned product does not contain melamine. However, the product is not tested for melamine.

This document has been created digitally and is valid without a signature.

BASF'S DISCLAIMERS AND WAIVERS ARE SET FORTH ON THE LAST PAGE

If you have any further questions or need additional support, please contact your BASF sales representative.

Kollidon® SR

PRD-No.: 30071321

Statement: Melamine free

This document, or any information provided herein does not constitute a legally binding obligation of BASF and has been prepared in good faith and is believed to be accurate as of the date of issuance. Unless expressly agreed otherwise in writing in a supply contract or other written agreement between you and BASF:

(a) To the fullest extent not prohibited by the applicable laws, BASF EXPRESSLY DISCLAIMS ALL OTHER REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES OF ANY KIND, WHETHER EXPRESS OR IMPLIED, WRITTEN OR ORAL, BY FACT OR LAW, INCLUDING ANY IMPLIED WARRANTIES, REPRESENTATIONS OR CONDITIONS OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, SATISFACTORY QUALITY, NON-INFRINGEMENT, AND ANY REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES, ARISING FROM STATUTE, COURSE OF DEALING OR USAGE OF TRADE and BASF HEREBY EXPRESSLY EXCLUDES AND DISCLAIMS ANY LIABILITY RESULTING FROM OR IN CONNECTION WITH THIS DOCUMENT OR ANY INFORMATION PROVIDED HEREIN, including, without limitation, any liability for any direct, consequential, special, or punitive damages relating to or arising therefrom, except in cases of:

- (i) death or personal injury to the extent caused by BASF's sole negligence;
- (ii) BASF's willful misconduct, fraud or fraudulent misrepresentation; or
- (iii) any matter in respect of which it would be unlawful for BASF to exclude or restrict liability under the applicable laws;

(b) Any information provided herein can be changed at BASF's sole discretion anytime and neither this document nor the information provided herein may be relied upon to satisfy from any and all obligations you may have to undertake your own inspections and evaluations;

(c) BASF rejects any obligation to, and will not, automatically update this document and any information provided herein, unless required by applicable law; and

(d) The user is responsible for confirming that the user has retrieved the most current version of this document from BASF as appropriate

Kollidon® SR
PRD-No.: 30071321

Statement: Phthalate

According to our knowledge of the production process, raw materials and equipment used, the above mentioned product does not contain phthalates. However, the product is not tested for phthalates.

This document has been created digitally and is valid without a signature.

BASF'S DISCLAIMERS AND WAIVERS ARE SET FORTH ON THE LAST PAGE

If you have any further questions or need additional support, please contact your BASF sales representative.

Kollidon® SR

PRD-No.: 30071321

Statement: Phthalate

This document, or any information provided herein does not constitute a legally binding obligation of BASF and has been prepared in good faith and is believed to be accurate as of the date of issuance. Unless expressly agreed otherwise in writing in a supply contract or other written agreement between you and BASF:

(a) To the fullest extent not prohibited by the applicable laws, BASF EXPRESSLY DISCLAIMS ALL OTHER REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES OF ANY KIND, WHETHER EXPRESS OR IMPLIED, WRITTEN OR ORAL, BY FACT OR LAW, INCLUDING ANY IMPLIED WARRANTIES, REPRESENTATIONS OR CONDITIONS OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, SATISFACTORY QUALITY, NON-INFRINGEMENT, AND ANY REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES, ARISING FROM STATUTE, COURSE OF DEALING OR USAGE OF TRADE and BASF HEREBY EXPRESSLY EXCLUDES AND DISCLAIMS ANY LIABILITY RESULTING FROM OR IN CONNECTION WITH THIS DOCUMENT OR ANY INFORMATION PROVIDED HEREIN, including, without limitation, any liability for any direct, consequential, special, or punitive damages relating to or arising therefrom, except in cases of:

- (i) death or personal injury to the extent caused by BASF's sole negligence;
- (ii) BASF's willful misconduct, fraud or fraudulent misrepresentation; or
- (iii) any matter in respect of which it would be unlawful for BASF to exclude or restrict liability under the applicable laws;

(b) Any information provided herein can be changed at BASF's sole discretion anytime and neither this document nor the information provided herein may be relied upon to satisfy from any and all obligations you may have to undertake your own inspections and evaluations;

(c) BASF rejects any obligation to, and will not, automatically update this document and any information provided herein, unless required by applicable law; and

(d) The user is responsible for confirming that the user has retrieved the most current version of this document from BASF as appropriate

Addressed to:

PFIZER INC



445 EASTERN POINT RD
06340-5157 GROTON
USA

Supply Chain Traceability Description

for Deliveries from Europe to North America


Kollidon® SR 20KG Plastic drums

PRD-No.: 30071321 | Article-No.: 52622254

The product mentioned above is manufactured at BASF SE, Carl-Bosch-Str. 38, 67056 Ludwigshafen, Germany including co-manufacturing at Schwarzwaldmilch GmbH, Bunsenstr. 2, 77652 Offenburg, Germany and tested for release at BASF SE, Carl Bosch Strasse 38, Ludwigshafen. From there the material is transported via truck to the site-contracted warehouse BASF SE, 67056 Warehouse-Verbund Ludwigshafen and Rhein-Main-Area, Germany (incl. TST GmbH, 67547 Worms/67550 Worms; Wetlog GmbH, 68169 Mannheim). Further, the product is transported by truck/rail and sea freight through the ports of Antwerp and Houston, then by truck/rail/barge to our regional distribution center NEDC Warehouse and Distribution Center, Edison, NJ, USA; and/or: Radiant Global Logistics, Bolton, ON, Canada. From there the product is delivered by truck/rail directly to consignees within North America.

There is no further analytical testing, repackaging and relabeling functions that the product undergoes before it is delivered to the consignee.

The following simplified process flow shows the path of the product mentioned above leaving BASF SE, Carl-Bosch-Str. 38, 67056 Ludwigshafen, Germany until it reaches consignees within North America. Please note that this scheme displays the planned general procedure. Deviations from this scheme might occur in individual cases as industry supply chains can be impacted by global, regional or local events. Also, in case you order through a distributor, there might be additional steps or warehouses before final delivery.



This is a computer generated document. It is valid without signature. This document including all the attachments is not intended for regulatory submissions. There is no automatic updating service.

BASF'S DISCLAIMERS AND WAIVERS ARE SET FORTH ON THE LAST PAGE

If you have any further questions or need additional support, please contact your BASF sales representative.

This document, or any information provided herein does not constitute a legally binding obligation of BASF and has been prepared in good faith and is believed to be accurate as of the date of issuance. Unless expressly agreed otherwise in writing in a supply contract or other written agreement between you and BASF:

(a) To the fullest extent not prohibited by the applicable laws, BASF EXPRESSLY DISCLAIMS ALL OTHER REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES OF ANY KIND, WHETHER EXPRESS OR IMPLIED, WRITTEN OR ORAL, BY FACT OR LAW, INCLUDING ANY IMPLIED WARRANTIES, REPRESENTATIONS OR CONDITIONS OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, SATISFACTORY QUALITY, NON-INFRINGEMENT, AND ANY REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES, ARISING FROM STATUTE, COURSE OF DEALING OR USAGE OF TRADE and BASF HEREBY EXPRESSLY EXCLUDES AND DISCLAIMS ANY LIABILITY RESULTING FROM OR IN CONNECTION WITH THIS DOCUMENT OR ANY INFORMATION PROVIDED HEREIN, including, without limitation, any liability for any direct, consequential, special, or punitive damages relating to or arising therefrom, except in cases of:

- (i) death or personal injury to the extent caused by BASF's sole negligence;
- (ii) BASF's willful misconduct, fraud or fraudulent misrepresentation; or
- (iii) any matter in respect of which it would be unlawful for BASF to exclude or restrict liability under the applicable laws;

(b) Any information provided herein can be changed at BASF's sole discretion anytime and neither this document nor the information provided herein may be relied upon to satisfy from any and all obligations you may have to undertake your own inspections and evaluations;

(c) BASF rejects any obligation to, and will not, automatically update this document and any information provided herein, unless required by applicable law; and

(d) The user is responsible for confirming that the user has retrieved the most current version of this document from BASF as appropriate

Kollidon® SR

PRD-No.: 30071321

Production Information

RegXcellence®

Production Information

Country of origin: Germany

Manufacturer

☒ BASF

Site (company name, location):

BASF SE, Carl-Bosch-Str. 38, 67056 Ludwigshafen, Germany

☒ Co-manufacturer (for certain process steps)

Site (company name, location):

Schwarzwaldmilch GmbH, Okenstr. 81a, 77652 Offenburg, Germany

☐ The product is manufactured by a contract manufacturer / 3rd party manufacturer according to the BASF specifications and requirements.

Original manufacturer:

Manufacturing process

☒ synthetic

☐ semi-synthetic following fermentation

☐ purification of natural products

☒ Others: co-processed excipient

Kollidon® SR

PRD-No.: 30071321

Production Information

RegXcellence®

This document, or any information provided herein does not constitute a legally binding obligation of BASF and has been prepared in good faith and is believed to be accurate as of the date of issuance. Unless expressly agreed otherwise in writing in a supply contract or other written agreement between you and BASF:

(a) To the fullest extent not prohibited by the applicable laws, BASF EXPRESSLY DISCLAIMS ALL OTHER REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES OF ANY KIND, WHETHER EXPRESS OR IMPLIED, WRITTEN OR ORAL, BY FACT OR LAW, INCLUDING ANY IMPLIED WARRANTIES, REPRESENTATIONS OR CONDITIONS OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, SATISFACTORY QUALITY, NON-INFRINGEMENT, AND ANY REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES, ARISING FROM STATUTE, COURSE OF DEALING OR USAGE OF TRADE and BASF HEREBY EXPRESSLY EXCLUDES AND DISCLAIMS ANY LIABILITY RESULTING FROM OR IN CONNECTION WITH THIS DOCUMENT OR ANY INFORMATION PROVIDED HEREIN, including, without limitation, any liability for any direct, consequential, special, or punitive damages relating to or arising therefrom, except in cases of:

- (i) death or personal injury to the extent caused by BASF's sole negligence;
- (ii) BASF's willful misconduct, fraud or fraudulent misrepresentation; or
- (iii) any matter in respect of which it would be unlawful for BASF to exclude or restrict liability under the applicable laws;

(b) Any information provided herein can be changed at BASF's sole discretion anytime and neither this document nor the information provided herein may be relied upon to satisfy from any and all obligations you may have to undertake your own inspections and evaluations;

(c) BASF rejects any obligation to, and will not, automatically update this document and any information provided herein, unless required by applicable law; and

(d) The user is responsible for confirming that the user has retrieved the most current version of this document from BASF as appropriate.

Technical Information




Kollidon[®] SR

Polyvinyl acetate and povidone based matrix sustained release excipient

August 2024 | Supersedes issue dated September 2015 | Last change DAWF-2024-0703

03_030728e-10/Page 1 of 9

® = Registered trademark of BASF in many countries.



 **BASF**

We create chemistry

Contents


	Page
1. Introduction	3
1.1 General	3
1.2 Chemical structure	3
1.3 Trivial name	3
2. Composition	3
3. Specifications and methods	3
3.1 Specification	3
3.2 IR-Spectra	3
4. Properties	4
5. Registration	5
5.1 Regulatory status	5
5.2 Drug Master File	5
5.3 Use of polyvinyl acetate in drugs and food	5
6. Applications	5
6.1 General Information	5
6.2 Propranolol substaisted release matrix tablets	7
6.3 Diclofenac substaisted release matrix tablets	7
6.4 Theophylline substaisted release matrix tablets	8
7. Stability & Safety	9
8. Available Articles and Packaging	9
9. Documents, Quality & Virtual Regulatory Information	9

1. Introduction

1.1 General

Kollidon® SR is a polyvinyl acetate and povidone based matrix retarding agent. It is particularly suitable for the manufacture of pH-independent sustained-release matrix tablets by direct compression. Polyvinyl acetate is a very plastic material that forms a coherent matrix even under low compression forces. When the tablets are introduced into gastric or intestinal fluid, the water soluble povidone is leached out to form pores allowing the active ingredient to diffuse. Kollidon® SR does not contain ionic groups and therefore does not show interactions with drug substances. The sustained-release properties are unaffected by ions or salts.

1.2 Chemical structure



1.3 Trivial name

Polyvinyl acetate/polyvinylpyrrolidone

2. Compositions

Kollidon® SR is a co-processed excipient consisting of polyvinyl acetate (from Kollicoat® SR 30 D) as a major component and povidone (from Kollidon® 30) as pore former. For the exact composition please refer to the product specification.

3. Specifications and methods

3.1 Specification

The current version of the product specification is available on [RegXcellence](#).

3.2 IR-spectra

The IR-spectra is measured in potassium bromide and a typical spectra is given in the following figure 1.




Fig. 1: IR-spectra of Kollidon® SR

4. Properties

Description	White or slightly yellowish, free-flowing powder.
Solubility	Insoluble in water (The povidone part is soluble but the polyvinyl acetate part is insoluble). It is very soluble in N-methylpyrrolidone.
Particle size distribution	The average particle size is about 50-70 µm.
Glass transition temperature	The glass transition temperature T _g of the anhydrous material is about 35 °C.
Bulk density	About 0.45 g/ml.
Flowability	Kollidon® SR has outstanding flow properties with a angle of repose well below 30 °. It can enhance the flowability of other components added for a tablet formulation.
Hygroscopicity	The water uptake is much less than that of povidone or copovidone. Figure 2 shows the water sorption and desorption isotherms at room temperature.




Fig. 2: Sorption isotherms of Kollidon® SR

Compressibility	Kollidon® SR has excellent compressibility and endows tablets with enormous hardness and low friability. This is due to the combination of the very plastic polyvinyl acetate and the also strongly binding povidone.
-----------------	---




Fig. 3: Hardness-compression force profile of propranolol sustained release tablets containing 50% of Kollidon® SR (2 Lots, Formulation see chapter 6.2)

5. Registration

5.1 Regulatory status

The product is a co-processed excipient. Thus monographs for Kollidon® SR do not exist.

5.2 Drug Master File

For registration purposes a US-DMF was prepared.

5.3 Use of polyvinyl acetate in drugs and food

Polyvinyl acetate is used in a variety of drugs for oral administration in numerous countries including Germany, France, Japan and USA. Polyvinyl acetate also is allowed in the food industry in several countries like Germany, USA and Japan.

6. Applications

6.1 General Information

Kollidon® SR can be used for the production of the following sustained release matrix dosage forms: Tablets, pellets and granules.

Different technologies to obtain such dosage forms can be applied: Direct compression, roller compaction, wet granulation and extrusion.

The excellent flowability and compressibility of Kollidon® SR makes this excipient particularly suitable for the manufacture of sustained release tablets obtained by **direct compression**.

The required content of Kollidon® SR in the tablet depends on the solubility of the active ingredient. The following table gives an information about the usual amounts of Kollidon® SR to obtain a sustained release during 12 – 24 hours.

Solubility of the active ingredient	Kollidon® SR in the tablet
Very slightly soluble to practically insoluble	15 – 25%
Sparingly soluble to slightly soluble	25 – 40%
Soluble to freely soluble	40 – 55%

The sustained release characteristics can be modified by varying the Kollidon® SR content in the formulation. Figure 4 shows the influence of the amount of Kollidon® SR on the release of caffeine as a example of a soluble active ingredient.




Fig. 4: Influence of the amount of Kollidon® SR on the drug release in a caffeine sustained release tablet (160 mg Caffeine)

In the case of slightly soluble or practically insoluble drug substances the release can be accelerated not only by reducing the content of Kollidon® SR but also by the addition of hydrophilic substances like lactose, Kollidon® 30 or Kollidon® CL-M which act as pore former.

Interesting and important properties of sustained release matrix tablets based on Kollidon® SR are the following:

1. The drug release is independent of the pH (see figure 5).
2. The drug release is independent of the ionic strength of the dissolution medium (see figure 5, addition of 2.5% of NaCl).
3. The drug release is independent of the usual compression force and tablet hardness (see figure 6).




Fig. 5: Influence of the pH and the ionic strength of the dissolution medium on the release of caffeine tablets (Caffeine + Kollidon® SR 1+1)

It is recommended to store the matrix tablets containing Kollidon® SR at temperatures below 30 °C and in tightly closed containers to avoid the uptake of humidity which could modify the release profile of formulations containing a higher percentage of Kollidon® SR.

In the following chapters three typical examples of soluble and practically insoluble active ingredients are given in form of sustained release tablets.

6.2

Propranolol Sustained Release Matrix Tablets

Formulation

	Parts by weight [g]	Composition [%]
Propranolol-HCl	160.0	49.23
Kollidon® SR	160.0	49.23
Silicon dioxide, colloidal	3.4	1.05
Magnesium stearate	1.6	0.49
Total	325.0	100.00

Manufacture

Tablet properties




Fig. 6: Propranolol sustained release tablets: Influence of the compression force on the drug release

6.3 Diclofenac Sustained Release Matrix Tablets

Formulation

	Weight	Percent
Diclofenac sodium	100 g	48.4
Kollidon® SR	100 g	48.4
Aerosil 200	3.4 g	1.6
Magnesium stearate	3.4 g	1.6

Manufacture

All ingredients are mixed, passed through a 0.8 mm sieve and pressed with a medium compression force on a rotary press.

Tablet properties

Diameter	8 mm
Weight	206 mg
Compression force	medium
Hardnes	195 N
Friability	<0.1%
Drug release	See Figure 7




Fig. 7: Dissolution of Diclofenac sustained release tablets

6.4

Theophylline Sustained Release Matrix Tablets

Formulation

	Parts by weight [g]	Composition [%]
Theophylline gran.	500.0	53.9
Kollidon® SR	200.0	21.6
Ludipress® LCE	225.0	24.2
Magnesium stearate	3.0	0.3
Total	928.0	100.00

Manufacture

All ingredients were passed through a 0.8 mm sieve, blended for 10 min in a Turbula mixer and then pressed on a rotary press.

Tablet properties

Diameter	19.0 x 8.5 mm (football shape)
Weight	928 mg
Compression force	11 kN
Hardness	172 N
Friability	<0.1%
Drug release	See Figure 8




Fig. 8: Dissolution of theophylline sustained release tablets

7. Stability & Safety

The product is typically stable for 48 months after date of production provided storage under recommended conditions. The actual retest period and storage conditions can be found in the document “Quality & Regulatory Product Summary” in [RegXcellence](#).

The actual version of the safety data sheet is accessible via [MyProductWorld](#) and send with every consignment.

8. Available Articles and Packaging

Product name	PRD number*	Article number	Packaging
Kollidon® SR	30071321	52622254	20 kg plastic drums

*BASF's commercial product number

9. Documents, Quality & Regulatory Information

Virtual Pharma Assistants



Meet your Virtual Pharma Assistants!
MyProductWorld, RegXcellence® and ZoomLab™

www.virtualpharmaassistants.com

[MyProductWorld](#): article numbers, sample order, safety data sheet, sustainability information

[RegXcellence](#): specification, compliance documents, regulatory product summary

[ZoomLab](#) formulation assistance to predict starting formulations and expedite drug formulations.

Disclaimer

This document and any information provided herein does not constitute a legally binding obligation of BASF and has been prepared in good faith and is believed to be accurate as of the date of issuance. Unless expressly agreed otherwise in writing in a supply contract or other written agreement between you and BASF:

- To the fullest extent not prohibited by the applicable laws, BASF EXPRESSLY DISCLAIMS ALL REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES OF ANY KIND, WHETHER EXPRESS OR IMPLIED, WRITTEN OR ORAL, BY FACT OR LAW, INCLUDING ANY IMPLIED WARRANTIES, REPRESENTATIONS OR CONDITIONS OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, SATISFACTORY QUALITY, NON-INFRINGEMENT, AND ANY REPRESENTATIONS, WARRANTIES, CONDITIONS OR GUARANTEES, ARISING FROM STATUTE, COURSE OF DEALING OR USAGE OF TRADE AND BASF HEREBY EXPRESSLY EXCLUDES AND DISCLAIMS ANY LIABILITY RESULTING FROM OR IN CONNECTION WITH THIS DOCUMENT OR ANY INFORMATION PROVIDED HEREIN, including, without limitation, any liability for any direct, consequential, special, or punitive damages relating to or arising therefrom, except in cases of (i) death or personal injury, (ii) BASF's or its agents and assistants willful misconduct, fraud or fraudulent misrepresentation or (iii) any matter in respect of which it would be unlawful for BASF to exclude or restrict liability under the applicable laws;
- Any information provided herein can be changed at BASF's sole discretion anytime and neither this document nor the information provided herein may be relied upon to satisfy any obligations you may have to undertake your own inspections and evaluation;
- BASF rejects any obligation to, and will not, automatically update this document and any information provided herein, unless required by applicable law; and
- You are responsible for confirming that you have retrieved the most current version of this document from BASF.

August 2024

Safety Data Sheet

Kollidon® SR

Revision date : 2025/08/28

Version: 6.0

Page: 1/10

(30071321/SDS_GEN_US/EN)

1. Identification

Product identifier used on the label

Kollidon® SR

Recommended use of the chemical and restriction on use

Recommended use*: pharmaceutical excipient

Unsuitable for use: Not intended for sale to or use by the general public.

* The "Recommended use" identified for this product is provided solely to comply with a Federal requirement and is not part of the seller's published specification. The terms of this Safety Data Sheet (SDS) do not create or infer any warranty, express or implied, including by incorporation into or reference in the seller's sales agreement.

Details of the supplier of the safety data sheet

Company:

BASF CORPORATION

100 Park Avenue

Florham Park, NJ 07932, USA

Telephone: +1 973 245-6000

Emergency telephone number

24 Hour Emergency Response Information

CHEMTREC: 1-800-424-9300

BASF HOTLINE: 1-800-832-HELP (4357)

Other means of identification

Synonyms:

Polymer based on: polyvinylacetate, 2Pyrrolidinone, 1-ethenyl-, homopolymer.

2. Hazards Identification

According to Regulation 2024 OSHA Hazard Communication Standard; 29 CFR Part 1910.1200

Classification of the product

Combustible Dust

Combustible Dust (1)

Combustible Dust

Label elements

Signal Word:

Safety Data Sheet

Kollidon® SR

Revision date: 2025/08/28
Version: 6.0

Page: 2/10
(30071321/SDS_GEN_US/EN)

Warning

Hazard Statement:

May form combustible dust concentration in air.

Hazards not otherwise classified

The product is under certain conditions capable of dust explosion.

3. Composition / Information on Ingredients

According to Regulation 2024 OSHA Hazard Communication Standard; 29 CFR Part 1910.1200

Under the referenced regulation, this product does not contain any components classified for health hazards above the relevant cut off value.

4. First-Aid Measures

Description of first aid measures

General advice:

Remove contaminated clothing.

If inhaled:

Keep patient calm, remove to fresh air.

If on skin:

Wash thoroughly with soap and water

If in eyes:

Wash affected eyes for at least 15 minutes under running water with eyelids held open.

If swallowed:

Rinse mouth and then drink 200-300 ml of water.

Most important symptoms and effects, both acute and delayed

Symptoms: No data available.

Indication of any immediate medical attention and special treatment needed

Note to physician

Treatment: Symptomatic treatment (decontamination, vital functions).

5. Fire-Fighting Measures

Extinguishing media

Safety Data Sheet

Kollidon® SR

Revision date: 2025/08/28
Version: 6.0

Page: 3/10
(30071321/SDS_GEN_US/EN)

Suitable extinguishing media:
foam, dry powder, water spray

Unsuitable extinguishing media for safety reasons:
water jet

Additional information:
Avoid whirling up the material/product because of the danger of dust explosion.

Special hazards arising from the substance or mixture

Hazards during fire-fighting:
carbon oxides, cyanides, nitrogen oxides
The substances/groups of substances mentioned can be released in case of fire.

Advice for fire-fighters

Protective equipment for fire-fighting:
Firefighters should be equipped with self-contained breathing apparatus and turn-out gear.

Wear a self-contained breathing apparatus.

Further information:

Dispose of fire debris and contaminated extinguishing water in accordance with official regulations.
Cool endangered containers with water-spray.

Dusty conditions may ignite explosively in the presence of an ignition source causing flash fire.

6. Accidental release measures

Further accidental release measures:

Avoid the formation and build-up of dust - danger of dust explosion. Dust in sufficient concentration can result in an explosive mixture in air. Handle to minimize dusting and eliminate open flame and other sources of ignition. Avoid dispersal of dust in the air (e.g. by clearing dusty surfaces with compressed air).

Personal precautions, protective equipment and emergency procedures

Use personal protective clothing. Information regarding personal protective measures, see section 8.
Avoid dust formation. Avoid contact with the skin, eyes and clothing.

Environmental precautions

Do not discharge into drains/surface waters/groundwater.

Methods and material for containment and cleaning up

For small amounts: Contain with dust binding material and dispose of.
For large amounts: Sweep/shovel up.
Dispose of absorbed material in accordance with regulations. Avoid raising dust.

Nonsparking tools should be used.

7. Handling and Storage

Precautions for safe handling

Safety Data Sheet

Kollidon® SR

Revision date: 2025/08/28
Version: 6.0

Page: 4/10
(30071321/SDS_GEN_US/EN)

Avoid dust formation. Provide exhaust ventilation if dust is formed.

Protection against fire and explosion:

Avoid dust formation. Dust in sufficient concentration can result in an explosive mixture in air. Handle to minimize dusting and eliminate open flame and other sources of ignition. Routine housekeeping should be instituted to ensure that dusts do not accumulate on surfaces. Dry powders can build static electricity charges when subjected to the friction of transfer and mixing operations. Provide adequate precautions, such as electrical grounding and bonding, or inert atmospheres. Refer to NFPA 660 (2025) Standard for Combustible Dust and Particulate Solids. NFPA 660 is a combination of Standards NFPA 61 (Agriculture and Food), NFPA 484 (Metals), NFPA 652 (Fundamentals of Combustible Dusts), NFPA 654 (Standard for the Prevention of Fire and Dust Explosions from the Manufacturing, Processing, and Handling of Combustible Particulate Solids), NFPA 65 (Sulfur), and NFPA 664 (Woodworking/Processing). Consult NFPA 660 standard for relevant commodity-specific and general safety information.

Conditions for safe storage, including any incompatibilities

No applicable information available.

Further information on storage conditions: Keep container tightly closed and dry.

8. Exposure Controls/Personal Protection

No substance specific occupational exposure limits known.

Advice on system design:

Provide local exhaust ventilation to control dusts/mists.

It is recommended that all dust control equipment such as local exhaust ventilation and material transport systems involved in handling of this product contain explosion relief vents or an explosion suppression system or an oxygen deficient environment. Ensure that dust-handling systems (such as exhaust ducts, dust collectors, vessels, and processing equipment) are designed in a manner to prevent the escape of dust into the work area (i.e., there is no leakage from the equipment). Use only appropriately classified electrical equipment and powered industrial trucks.

Personal protective equipment

Respiratory protection:

Breathing protection if dusts are formed. Wear a NIOSH-certified (or equivalent) particulate respirator.

Hand protection:

Chemical resistant protective gloves

Eye protection:

Wear safety goggles (chemical goggles) if there is potential for airborne dust exposures.

Body protection:

Body protection must be chosen based on level of activity and exposure.

General safety and hygiene measures:

Handle in accordance with good industrial hygiene and safety practice. Wearing of closed work clothing is recommended. Wear protective clothing as necessary to minimize contact. No eating, drinking, smoking or tobacco use at the place of work. Hands and/or face should be washed before breaks and at the end of the shift. Store work clothing separately.

Safety Data Sheet

Kollidon® SR

Revision date: 2025/08/28
Version: 6.0

Page: 5/10
(30071321/SDS_GEN_US/EN)

9. Physical and Chemical Properties

Physical state:	solid
Form:	powder
Odour:	almost odourless
Odour threshold:	not determined
Colour:	white to slightly yellow
pH value:	3.5 - 5.5 (100 g/l, 20 °C)
melting range:	>= 130 °C The substance / product decomposes.
Freezing point:	No data available.
Boiling point:	not applicable
Sublimation point:	No applicable information available.
Flash point:	not applicable, the product is a solid
Flammability:	not highly flammable (VDI 2263, sheet 1, 1.1 (May 1990))
Lower explosion limit:	For solids not relevant for classification and labelling.
Upper explosion limit:	For solids not relevant for classification and labelling.
SADT:	Not a substance liable to self-decomposition according to UN transport regulations, class 4.1.
Vapour pressure:	negligible
Density:	No information is available for the absolute density. Instead the bulk density was determined as a more relevant value.
Bulk density:	400 kg/m ³
Relative vapour density:	The product is a non-volatile solid.
Partitioning coefficient n-octanol/water (log Pow):	not applicable
Thermal decomposition:	305 °C (DSC (DIN 51007)) It is not a self-decompositionable substance.
Viscosity, dynamic:	No data available.
Viscosity, kinematic:	not applicable, the product is a solid
Solubility in water:	dispersible, of low solubility
Solubility (quantitative):	No applicable information available.
Solubility (qualitative):	readily soluble solvent(s): N-Methylpyrrolidone
Molecular weight:	No data available.
Evaporation rate:	The product is a non-volatile solid.

Particle characteristics

Particle size distribution: typically > 60 µm (D50, Volumetric Distribution, ISO 13320-1)

10. Stability and Reactivity

Reactivity

No hazardous reactions if stored and handled as prescribed/indicated.

Corrosion to metals:

Safety Data Sheet

Kollidon® SR

Revision date: 2025/08/28
Version: 6.0

Page: 6/10
(30071321/SDS_GEN_US/EN)

Corrosive effects to metal are not anticipated.

Oxidizing properties:

Based on its structural properties the product is not classified as oxidizing.

Chemical stability

The product is stable if stored and handled as prescribed/indicated.

Possibility of hazardous reactions

Dust can form an explosive mixture with air.

Conditions to avoid

Avoid dust formation. Avoid electro-static charge. Avoid all sources of ignition: heat, sparks, open flame.

Incompatible materials

strong alkalies

Hazardous decomposition products

Decomposition products:

Hazardous decomposition products: No hazardous decomposition products if stored and handled as prescribed/indicated.

Thermal decomposition:

305 °C (DSC (DIN 51007))

It is not a self-decompositionable substance.

11. Toxicological information

Primary routes of exposure

Routes of entry for solids and liquids are ingestion and inhalation, but may include eye or skin contact. Routes of entry for gases include inhalation and eye contact. Skin contact may be a route of entry for liquefied gases.

Acute Toxicity/Effects

Acute toxicity

Assessment of acute toxicity: Virtually nontoxic after a single ingestion. The product has not been tested. The statement has been derived from substances/products of a similar structure or composition.

Oral

Type of value: LD50

Species: rat

Value: > 2,000 mg/kg (OECD Guideline 423)

Inhalation

No data available.

Type of value: ATE

Value: > 5 mg/l

Determined for dust

Safety Data Sheet

Kollidon® SR

Revision date: 2025/08/28
Version: 6.0

Page: 7/10
(30071321/SDS_GEN_US/EN)

Dermal

No data available.

Type of value: ATE

Value: > 5,000 mg/kg

Assessment other acute effects

Assessment of STOT single:

Based on available data, the classification criteria are not met.

Irritation / corrosion

Assessment of irritating effects: Not irritating to the skin. Not irritating to the eyes. The product has not been tested. The statement has been derived from substances/products of a similar structure or composition.

Skin

Species: rabbit

Result: non-irritant

Method: OECD Guideline 404

Eye

Species: rabbit

Result: non-irritant

Method: OECD Guideline 405

Sensitization

Assessment of sensitization: No data available.

Aspiration Hazard

not applicable

Chronic Toxicity/Effects

Repeated dose toxicity

Assessment of repeated dose toxicity: No data available.

Genetic toxicity

Assessment of mutagenicity: The substance was not mutagenic in bacteria. The substance was not mutagenic in mammalian cell culture. The substance was not mutagenic in a test with mammals. The product has not been tested. The statement has been derived from substances/products of a similar structure or composition.

Carcinogenicity

Assessment of carcinogenicity: No data available concerning carcinogenic effects.

Reproductive toxicity

Assessment of reproduction toxicity: Repeated oral uptake of the substance did not cause damage to the reproductive organs. The product has not been tested. The statement has been derived from substances/products of a similar structure or composition.

Teratogenicity

Assessment of teratogenicity: No indications of a developmental toxic / teratogenic effect were seen in animal studies. The product has not been tested. The statement has been derived from substances/products of a similar structure or composition.

Safety Data Sheet

Kollidon® SR

Revision date: 2025/08/28
Version: 6.0

Page: 8/10
(30071321/SDS_GEN_US/EN)

Other Information

The product has not been tested. The statements on toxicology have been derived from the properties of the individual components.

12. Ecological Information

Toxicity

Aquatic toxicity

Assessment of aquatic toxicity:

There is a high probability that the product is not acutely harmful to aquatic organisms. The inhibition of the degradation activity of activated sludge is not anticipated when introduced to biological treatment plants in appropriate low concentrations. The product has not been tested. The statement has been derived from substances/products of a similar structure or composition.

Aquatic invertebrates

EC50 (48 h) > 100 mg/l, Daphnia magna (OECD Guideline 202, part 1, static)

The product has low solubility in the test medium. An eluate has been tested.

Assessment of terrestrial toxicity

No data available.

Microorganisms/Effect on activated sludge

Toxicity to microorganisms

OECD Guideline 209 aerobic

activated sludge, domestic, adapted/EC20 (30 min): > 1,000 mg/l

Persistence and degradability

Assessment biodegradation and elimination (H2O)

Poorly biodegradable. The product has not been tested. The statement has been derived from substances/products of a similar structure or composition.

Elimination information

0 - 10 % BOD of the ThOD (28 d) (OECD Guideline 301 F) (aerobic, activated sludge, domestic)

Bioaccumulative potential

Assessment bioaccumulation potential

Significant accumulation in organisms is not to be expected.

Mobility in soil

Assessment transport between environmental compartments

No data available.

No data available.

Additional information

Add. remarks environm. fate & pathway:

The product has not been tested. The statements on environmental fate and pathway have been derived from the properties of the individual components.

Safety Data Sheet

Kollidon® SR

Revision date: 2025/08/28
Version: 6.0

Page: 9/10
(30071321/SDS_GEN_US/EN)

Other ecotoxicological advice:

The product has not been tested. The statements on ecotoxicology have been derived from the properties of the individual components.

13. Disposal considerations

Waste disposal of substance:

Do not discharge into drains/surface waters/groundwater. Dispose of in accordance with national, state and local regulations.

Container disposal:

Uncontaminated packaging can be re-used. Packs that cannot be cleaned should be disposed of in the same manner as the contents.

14. Transport Information

Land transport

USDOT

Not classified as a dangerous good under transport regulations

Sea transport

IMDG

Not classified as a dangerous good under transport regulations

Air transport

IATA/ICAO

Not classified as a dangerous good under transport regulations

15. Regulatory Information

Federal Regulations

Registration status:

Pharma TSCA, US released / exempt

Chemical TSCA, US released / listed

Chemical TSCA, US

All substances are TSCA listed and active.

EPCRA 311/312 (Hazard categories): Refer to SDS section 2 for GHS hazard classes applicable for this product.

NFPA Hazard codes:

Health: 0 Fire: 1 Reactivity: 0 Special:

HMIS III rating

Health: 0 Flammability: 1 Physical hazard: 0

Safety Data Sheet

Kollidon® SR

Revision date: 2025/08/28
Version: 6.0

Page: 10/10
(30071321/SDS_GEN_US/EN)

16. Other Information

SDS Prepared by:

BASF NA Product Regulations

SDS Prepared on: 2025/08/28

We support worldwide Responsible Care® initiatives. We value the health and safety of our employees, customers, suppliers and neighbors, and the protection of the environment. Our commitment to Responsible Care is integral to conducting our business and operating our facilities in a safe and environmentally responsible fashion, supporting our customers and suppliers in ensuring the safe and environmentally sound handling of our products, and minimizing the impact of our operations on society and the environment during production, storage, transport, use and disposal of our products.

Kollidon® SR is a registered trademark of BASF Corporation or BASF SE

IMPORTANT: WHILE THE DESCRIPTIONS, DESIGNS, DATA AND INFORMATION CONTAINED HEREIN ARE PRESENTED IN GOOD FAITH AND BELIEVED TO BE ACCURATE, IT IS PROVIDED FOR YOUR GUIDANCE ONLY. BECAUSE MANY FACTORS MAY AFFECT PROCESSING OR APPLICATION/USE, WE RECOMMEND THAT YOU MAKE TESTS TO DETERMINE THE SUITABILITY OF A PRODUCT FOR YOUR PARTICULAR PURPOSE PRIOR TO USE. NO WARRANTIES OF ANY KIND, EITHER EXPRESSED OR IMPLIED, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ARE MADE REGARDING PRODUCTS DESCRIBED OR DESIGNS, DATA OR INFORMATION SET FORTH, OR THAT THE PRODUCTS, DESIGNS, DATA OR INFORMATION MAY BE USED WITHOUT INFRINGING THE INTELLECTUAL PROPERTY RIGHTS OF OTHERS. IN NO CASE SHALL THE DESCRIPTIONS, INFORMATION, DATA OR DESIGNS PROVIDED BE CONSIDERED A PART OF OUR TERMS AND CONDITIONS OF SALE. FURTHER, YOU EXPRESSLY UNDERSTAND AND AGREE THAT THE DESCRIPTIONS, DESIGNS, DATA, AND INFORMATION FURNISHED BY OUR COMPANY HEREUNDER ARE GIVEN GRATIS AND WE ASSUME NO OBLIGATION OR LIABILITY FOR THE DESCRIPTION, DESIGNS, DATA AND INFORMATION GIVEN OR RESULTS OBTAINED, ALL SUCH BEING GIVEN AND ACCEPTED AT YOUR RISK. Any other intended applications should be discussed with the manufacturer. Corresponding occupational protection measurements must be followed.

Date / Revised: 2025/08/28
Date / Previous version: 2024/05/20

Version: 6.0
Previous version: 5.0

END OF DATA SHEET